



PRODUCT RANGE

A selection of the REAL-chairs



Valid from 080601

Last change 10 03 30

Article number 803382

A BRIEF DESCRIPTION OF THE COMPANY

Mercado Medic AB is engaged in the business of developing, selling and service of special chairs for all people with functional disabilities. We also produce specially adapted chairs to various clinics, offices and industries. Our head quarters are situated in Sweden and our export office, which distributes our products to the rest of Europe, in The Netherlands.



THE REAL-SYSTEM

The REAL-system is a Swedish quality system chair with a modular construction for simple individual adjustment.

With a variety of seats, back-rests, wheels, height range and accessories you can simply combine to ensure a functional and comfortable chair. It is suited perfectly to be used at home as well as in the office.





REAL 9000 PLUS

Sit in a REAL 9000 PLUS at home and at work. It is easy to adjust position for whatever you want to do. This model has a gas-regulated height adjustment of 20 cm. When the brakes are on, the chair stands firmly in position when you want to get up or transfer. The easy rolling wheels make mobility easy.

REAL 9100 PLUS EL

This model is equipped with electric height adjustment of an impressive 28 cm. By pressing the button you can raise or lower the automatic seat to the desired height. Brakes and the adjustable angle of the seat and backrest are also available with an electric regulation device.



REAL 9700/9800 PLUS COXIT

REAL 9000 PLUS is also available with a Coxit seat. The Coxit seat is divided with individually adjustable flaps. They are attached without any dust trapping recesses. This model is best suited for people with various types of orthopaedic injuries. For chairs with an adjustable seat angle the paddle levers are situated on at the back of the seat. On the Coxit chair the paddle lever is situated on the flap.



MADE IN SWEDEN

Manufacturing takes place in house, and all key Components are manufactured on modern tooling machines at head office on the island of Lidingö, Sweden.

1. Handmade upholstery
2. Comfortable armrests
3. Easy access handles
4. Safe and stable brakes
5. Stable frame
6. SKF castors

REAL[®] 6100



REAL 9300/9400 EL PLUS CHILD

Children with various disabilities are able to sit properly at a table With the REAL 9300 or 9400 EL PLUS. The adjustment functions has been developed in consultation with medical experts and provide active functional sitting posture. The REAL-chair is available with manual or electric lift as well as a number of functional options.



REAL 9200 TWIN/EL

is a module-chair intended as a means for moving about for somewhat heavier users. It is equipped with either double electric- or double gas spring actuator, a double back-mechanism and a reinforced X-base.

The electric model (REAL 9200 TWIN EL) fits Excellent as a “standing up aid” and has a capacity of up to 275 kg user weight. The manual model (REAL 9200 TWIN) has a capacity of up to 180 kg user weight.

REAL 6100 MOBILE

Is an electrically powered wheelchair designed for indoor use. It is constructed to provide the user maximum freedom of movement at home, in the office or at school. The range of electric height adjustment is an impressive 28 cm. This means that one has an extensive freedom of movement even in elevated positions. REAL 6100 is available with a number of useful options as well as a great variety of functional accessories.





REAL 9500/9600 EL

Are the chairs that are especially adapted for people with weak muscles who need to relieve their hips, knees or ankles while working in a standing position at home or at work. These chairs are equipped with, among other things, a H-frame, a saddle seat and an elevated seat position.

REAL 2000

Are included in a series of stand-supporting chairs to be used in standing-sitting working positions. By that we mean that you, while doing your work, have the option of standing with support. It is easier to avoid professional injuries caused by extensive static load by dividing the bodyweight between legs, seat area and back-support. These chairs can be ordered with a variety of saddle-seats, back-supports and seat levels as well as with or without wheels.



REAL 2624

Sitting in a wrong position during long working hours tires our bodies and may often cause wear damage. Weariness and tension in muscles are often the results of a prolonged sitting in a bent-forward position. That is why Mercado Medic has, in co-operation with physio- and occupational therapists, designed REAL 2621 and REAL 2624 where one can sit with a chest support instead of a back support.

MERCADO MEDIC UK

HENRIK SELLBERG

Managing Director

VISITING ADDRESS

Mercado Medic UK Ltd
Unit 7, Exis Court
Veasey Close
Nuneaton
CV11 6RT

CONTACT INFORMATION

Postal Address PO Box 4818, Warwick, CV34 9GB
Tel +44 (0) 844 561 6874
Fax +44 (0) 844 561 6875
Mobile +44 (0) 7855 794 757

Internet www.mercado-medic.co.uk
E-mail Sales@mercado-medic.co.uk
Enquiries info@mercado-medic.co.uk



Mercado Medic UK Ltd
PO Box 4818
Warwick, CV34 9GB

Phone & Fax
Ph +44-(0)844 561 6874
Fax +44-(0)844 561 6875

Email & Internet
Email info@mercado-medic.co.uk
Internet www.mercado-medic.co.uk